

Product Datasheet

ATP6V0A2 Antibody (orb1245642)

Description

ATP6V0A2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ATP6V0A2.

Target

ATP6V0A2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ATP6V0A2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ATP6V0A2 antibody can be used for detection of ATP6V0A2 by ELISA at 1:312500. ATP6V0A2 antibody can be used for detection of ATP6V0A2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

98 kDa

Uniprot ID

[Q9Y487](#)

NCBI


[NP_036595](#)

Dilution Range


ATP6V0A2 antibody can be used for detection of ATP6V0A2 by ELISA at 1:312500. ATP6V0A2 antibody can be used for detection of ATP6V0A2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on human kidney at ...



Antibody used in IHC on Human small inte...



Antibody used in IHC on Human kidney at ...