

Product Datasheet

SLC6A2 Antibody (orb1245675)

Description

SLC6A2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC6A2.

Target

SLC6A2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC6A2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC6A2 antibody can be used for detection of SLC6A2 by ELISA at 1:312500. SLC6A2 antibody can be used for detection of SLC6A2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

69 kDa

Uniprot ID

[P23975](#)

NCBI

[NP_001034](#)

Dilution Range

SLC6A2 antibody can be used for detection of SLC6A2 by ELISA at 1:312500. SLC6A2 antibody can be used for detection of SLC6A2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90kDa

40kDa

31kDa

22kDa

Antibody used in WB on Human Heart at 0...

SLC6A2

90kDa

65kDa

40kDa

29kDa

22kDa

Rabbit Anti SLC6A2
Sample Type Human Fetal Heart
Antibody Concentration 1µg/mL

Antibody used in WB on Hum. Fetal Heart ...

SLC6A2

90kDa

65kDa

40kDa

29kDa

22kDa

Rabbit Anti SLC6A2
Sample Type Human Fetal Heart
Antibody Concentration 1µg/mL

Antibody used in WB on Hum. Fetal Heart ...