


## Product Datasheet

### ALB Antibody (orb1245686)

## Description

ALB Antibody



## Species/Host

Rabbit

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

ELISA, WB

## Applications

Antibody used in IHC on Human Kidney.

## Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human Albumin.

## Target

ALB

## Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

## Form/Appearance

Liquid

## Concentration

batch dependent

## Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store Albumin antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

## Note

For research use only

## Application notes

Albumin antibody can be used for detection of Albumin by ELISA at 1:312500. Albumin antibody can be used for detection of Albumin by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

## Clonality

Polyclonal

## MW

67 kDa

## Uniprot ID

[P02768](#)

## NCBI

[NP\\_000468](#)

## Dilution Range

Albumin antibody can be used for detection of Albumin by ELISA at 1:312500. Albumin antibody can be used for detection of Albumin by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

## Expiration Date

12 months from date of receipt.