

Product Datasheet

RHOU Antibody (orb1245719)

Description

RHOU Antibody

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RHOU.
Target	RHOU
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RHOU antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	RHOU antibody can be used for detection of RHOU by ELISA at 1:1562500. RHOU antibody can be used for detection of RHOU by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	28 kDa
Uniprot ID	Q7L0Q8
NCBI	NP_067028
Dilution Range	RHOU antibody can be used for detection of RHOU by ELISA at 1:1562500. RHOU antibody can be used for detection of RHOU by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Muscle at 0...