

## Product Datasheet

### A1CF Antibody (orb1245736)

## Description

A1CF Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human A1CF.

Antibody used in WB on Human HepG2 at 0....

### Target

A1CF

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store A1CF antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

A1CF antibody can be used for detection of A1CF by ELISA at 1:1562500. A1CF antibody can be used for detection of A1CF by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

65 kDa

### Uniprot ID

[Q9NQ94](#)

### NCBI

[NP\\_620311](#)

### Dilution Range

A1CF antibody can be used for detection of A1CF by ELISA at 1:1562500. A1CF antibody can be used for detection of A1CF by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.

90 kDa

65 kDa

31 kDa

22 kDa

A1CF

90

65

40

31

22

Rabbit Anti-A1CF  
Sample Type: Human Adult  
Placenta  
Antibody Concentration: Liquid

Antibody used in WB on Hum. Adult Placen...