

Product Datasheet

TNRC6A Antibody (orb1245749)

Description

TNRC6A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
MCF-7 at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TNRC6A.

Target

TNRC6A

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TNRC6A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TNRC6A antibody can be used for detection of TNRC6A by ELISA at 1:62500. TNRC6A antibody can be used for detection of TNRC6A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

210 kDa

Uniprot ID

[Q8NDV7](#)

NCBI

[NP_055309](#)

Dilution Range

TNRC6A antibody can be used for detection of TNRC6A by ELISA at 1:62500. TNRC6A antibody can be used for detection of TNRC6A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.