

Product Datasheet

KLHDC8A Antibody (orb1245752)

Description

KLHDC8A Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KLHDC8A.
Target	KLHDC8A
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store KLHDC8A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	KLHDC8A antibody can be used for detection of KLHDC8A by ELISA at 1:1562500. KLHDC8A antibody can be used for detection of KLHDC8A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q8IYD2
NCBI	NP_060673
Dilution Range	KLHDC8A antibody can be used for detection of KLHDC8A by ELISA at 1:1562500. KLHDC8A antibody can be used for detection of KLHDC8A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

Antibody
used in WB
on Human
HepG2 at
0....