

Product Datasheet

CYP11A1 antibody (PE) (orb124576)

Description

Rabbit polyclonal to CYP11A1 (PE).

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

PE

Tested

FC

Applications
Immunogen

KLH conjugated synthetic peptide derived from human CYP11A1/P450SCC

Target

CYP11A1

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Application notes

orb124576 is the PE conjugated form of orb5936

Isotype

IgG

Clonality

Polyclonal

MW

53/57kDa

Uniprot ID
P05108
Dilution Range

Flow-Cyt=,IF=1:100-500

Expiration Date

12 months from date of receipt.