

Product Datasheet

CDC5L Antibody (orb1245801)

Description

CDC5L Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
721_B at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CDC5L.

Target

CDC5L

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CDC5L antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CDC5L antibody can be used for detection of CDC5L by ELISA at 1:312500. CDC5L antibody can be used for detection of CDC5L by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

92 kDa

Uniprot ID

[Q99459](#)

NCBI

[NP_001244](#)

Dilution Range

CDC5L antibody can be used for detection of CDC5L by ELISA at 1:312500. CDC5L antibody can be used for detection of CDC5L by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.