

Product Datasheet

TUBB2A Antibody (orb1245812)

Description

TUBB2A Antibody

Species/Host

Rabbit

Reactivity

Drosophila, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TUBB2A.

Antibody used in WB on Human HepG2 at 1

Target

TUBB2A

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TUBB2A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TUBB2A antibody can be used for detection of TUBB2A by ELISA at 1:312500. TUBB2A antibody can be used for detection of TUBB2A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

49 kDa


Uniprot ID
[Q13885](#)
NCBI
[NP_001060](#)
Dilution Range

TUBB2A antibody can be used for detection of TUBB2A by ELISA at 1:312500. TUBB2A antibody can be used for detection of TUBB2A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90 kDa
66 kDa
40 kDa
31 kDa
22 kDa



Antibody used in IHC on Small intestine,...