

Product Datasheet

FOXK2 Antibody (orb1245827)

Description

FOXK2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human FOXK2.

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human FOXK2.

Target

FOXK2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store FOXK2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

FOXK2 antibody can be used for detection of FOXK2 by ELISA at 1:1562500. FOXK2 antibody can be used for detection of FOXK2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

69 kDa

Uniprot ID

[Q01167](#)

NCBI

[NP_004505](#)

Dilution Range

FOXK2 antibody can be used for detection of FOXK2 by ELISA at 1:1562500. FOXK2 antibody can be used for detection of FOXK2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.


Expiration Date

12 months from date of receipt.




Antibody used in WB on Human 293T at 0.2...

FOXK2



Antibody used in WB on Hum. Fetal Brain ...

FOXK2



Antibody used in WB on Human 721_B at 1 ...