

Product Datasheet

SLC26A10 Antibody (orb1245847)

Description	SLC26A10 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of SLC26A10.	Antibody used in WB on Human HepG2 at 0....
Target	SLC26A10	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC26A10 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	SLC26A10 antibody can be used for detection of SLC26A10 by ELISA at 1:62500. SLC26A10 antibody can be used for detection of SLC26A10 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	60 kDa	
Uniprot ID	Q8NG04	
NCBI	NP_597996	
Dilution Range	SLC26A10 antibody can be used for detection of SLC26A10 by ELISA at 1:62500. SLC26A10 antibody can be used for detection of SLC26A10 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	