

Product Datasheet

CPXCR1 Antibody (orb1245853)

Description

CPXCR1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CPXCR1.

Target

CPXCR1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CPXCR1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CPXCR1 antibody can be used for detection of CPXCR1 by ELISA at 1:312500. CPXCR1 antibody can be used for detection of CPXCR1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

[Q8N123](#)

NCBI


[NP_149037](#)

Dilution Range


CPXCR1 antibody can be used for detection of CPXCR1 by ELISA at 1:312500. CPXCR1 antibody can be used for detection of CPXCR1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human thymus at ...



Antibody used in IHC on Human testis at ...



Antibody used in WB on Human Fetal Brain...