

Product Datasheet

ZNF185 Antibody (orb1245943)

Description

ZNF185 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF185.

Target

ZNF185

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF185 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZNF185 antibody can be used for detection of ZNF185 by ELISA at 1:12500. ZNF185 antibody can be used for detection of ZNF185 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

73 kDa

Uniprot ID

[A4FTV3](#)

NCBI


[NP_009081](#)

Dilution Range


ZNF185 antibody can be used for detection of ZNF185 by ELISA at 1:12500. ZNF185 antibody can be used for detection of ZNF185 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in WB on Human HeLa at 0.2...



Antibody used in WB on Hum. Fetal Brain ...



Antibody used in WB on Human MCF7 at 1 u...