

Product Datasheet

TBX18 Antibody (orb1246004)

Description

TBX18 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human TBX18.

Antibody used in WB on Human Jurkat 0.2-...

Target

TBX18

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TBX18 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in IHC on Human Breast at ...

Note

For research use only

Application notes

TBX18 antibody can be used for detection of TBX18 by ELISA at 1:312500. TBX18 antibody can be used for detection of TBX18 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

65 kDa

Uniprot ID

[O95935](#)

NCBI

[XP_943361](#)

Dilution Range

TBX18 antibody can be used for detection of TBX18 by ELISA at 1:312500. TBX18 antibody can be used for detection of TBX18 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

