

Product Datasheet

GJA9 Antibody (orb1246029)


Description

GJA9 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GJA9.
Target	GJA9
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GJA9 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GJA9 antibody can be used for detection of GJA9 by ELISA at 1:62500. GJA9 antibody can be used for detection of GJA9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	P57773
NCBI	NP_110399
Dilution Range	GJA9 antibody can be used for detection of GJA9 by ELISA at 1:62500. GJA9 antibody can be used for detection of GJA9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Brain at 0....



Antibody used in WB on Hum. Adult Placen...



Antibody used in WB on Hum. Fetal Liver ...