

## Product Datasheet

### ZNF75A Antibody (orb1246066)

## Description

ZNF75A Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

Antibody  
used in WB  
on Human  
MCF-7 at 0....

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF75A.

### Target

ZNF75A

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF75A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

ZNF75A antibody can be used for detection of ZNF75A by ELISA at 1:12500. ZNF75A antibody can be used for detection of ZNF75A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

35 kDa

### Uniprot ID

[Q96N20](#)

### NCBI

[NP\\_694573](#)

### Dilution Range

ZNF75A antibody can be used for detection of ZNF75A by ELISA at 1:12500. ZNF75A antibody can be used for detection of ZNF75A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.