

Product Datasheet

ZFP3 Antibody (orb1246067)

Description

ZFP3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HepG2 at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZFP3.

Target

ZFP3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZFP3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZFP3 antibody can be used for detection of ZFP3 by ELISA at 1:1562500. ZFP3 antibody can be used for detection of ZFP3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

58 kDa

Uniprot ID

[Q96NJ6](#)

NCBI

[NP_694563](#)

Dilution Range

ZFP3 antibody can be used for detection of ZFP3 by ELISA at 1:1562500. ZFP3 antibody can be used for detection of ZFP3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.