

Product Datasheet


SLC18A1 Antibody (orb1246128)

Description


SLC18A1 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC18A1.
Target	SLC18A1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC18A1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SLC18A1 antibody can be used for detection of SLC18A1 by ELISA at 1:62500. SLC18A1 antibody can be used for detection of SLC18A1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P54219
NCBI	NP_003044
Dilution Range	SLC18A1 antibody can be used for detection of SLC18A1 by ELISA at 1:62500. SLC18A1 antibody can be used for detection of SLC18A1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

Antibody used in WB on Human 293T at 0.2...



Antibody used in IHC on Human brain at 2...



Antibody used in IHC on Human brain.