

Product Datasheet

KCNN2 Antibody (orb1246140)


Description

KCNN2 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KCNN2.
Target	KCNN2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store KCNN2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	KCNN2 antibody can be used for detection of KCNN2 by ELISA at 1:62500. KCNN2 antibody can be used for detection of KCNN2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	26 kDa
Uniprot ID	Q6PJ10
NCBI	NP_740721
Dilution Range	KCNN2 antibody can be used for detection of KCNN2 by ELISA at 1:62500. KCNN2 antibody can be used for detection of KCNN2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

40kDa
29kDa
22kDa

Antibody used in WB on Human Muscle at 0...



Green: KCNN2

See HC-1 Data and Customer Feedback for more information

Antibody used in IHC on Rhesus macaque s...

KCNN2

65
40
29
22

Rabbit Anti-KCNN2
Sample Type: Human Fetal
Muscle
Antibody Concentration: 1µg/mL

Antibody used in WB on Hum. Fetal Muscle...