

Product Datasheet


KCNJ1 Antibody (orb1246141)

Description

KCNJ1 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KCNJ1.
Target	KCNJ1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store KCNJ1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	KCNJ1 antibody can be used for detection of KCNJ1 by ELISA at 1:312500. KCNJ1 antibody can be used for detection of KCNJ1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	P48048
NCBI	NP_000211
Dilution Range	KCNJ1 antibody can be used for detection of KCNJ1 by ELISA at 1:312500. KCNJ1 antibody can be used for detection of KCNJ1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

KCNJ1



See Immunoblot 2700 and Customer Feedback for more information

Antibody used in WB on Mouse renal epith...



Antibody used in WB on Human HepG2 at 0....