

Product Datasheet

KCTD7 Antibody (orb1246142)


Description

KCTD7 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KCTD7.
Target	KCTD7
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store KCTD7 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	KCTD7 antibody can be used for detection of KCTD7 by ELISA at 1:62500. KCTD7 antibody can be used for detection of KCTD7 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	Q96MP8
NCBI	NP_694578
Dilution Range	KCTD7 antibody can be used for detection of KCTD7 by ELISA at 1:62500. KCTD7 antibody can be used for detection of KCTD7 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Transfected 293T ...



Antibody used in WB on Hum. Adult Placen...



Antibody used in WB on Hum. Fetal Lung a...