

Product Datasheet

ErbB 2 antibody (PE) (orb124615)

Description

Rabbit polyclonal to ErbB 2 (PE).

Species/Host

Rabbit

Reactivity

Human

Conjugation

PE

Tested Applications

ICC, IF

Immunogen

KLH conjugated synthetic peptide derived from human HER2

Target

ERBB2

Preservatives

0.01M PBS(pH7.4)

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Application notes

orb124615 is the PE conjugated form of orb10372

Isotype

IgG

Clonality

Polyclonal

MW

185kDa

Uniprot ID
P04626
Dilution Range

IF=1:100-500,ICC=1:100

Expiration Date

12 months from date of receipt.