

Product Datasheet

TRPC4 Antibody (orb1246151)


Description

TRPC4 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TRPC4.
Target	TRPC4
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TRPC4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	TRPC4 antibody can be used for detection of TRPC4 by ELISA at 1:312500. TRPC4 antibody can be used for detection of TRPC4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	112 kDa
Uniprot ID	Q9UBN4
NCBI	NP_057263
Dilution Range	TRPC4 antibody can be used for detection of TRPC4 by ELISA at 1:312500. TRPC4 antibody can be used for detection of TRPC4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human 293T at 0.2...



Antibody used in IHC on Human Placenta a...



Antibody used in WB on Human MCF7 at 1 u...