

Product Datasheet

CNBP Antibody (orb1246172)

Description

CNBP Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CNBP.

Target

CNBP

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CNBP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CNBP antibody can be used for detection of CNBP by ELISA at 1:2500. CNBP antibody can be used for detection of CNBP by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

19 kDa

Uniprot ID

[P62633](#)

NCBI

[NP_003409](#)

Dilution Range


CNBP antibody can be used for detection of CNBP by ELISA at 1:2500. CNBP antibody can be used for detection of CNBP by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on human primary myo...



Antibody used in WB on Human brain at 0....