

Product Datasheet

PRDM8 Antibody (orb1246187)

Description

PRDM8 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human PRDM8.

Target

PRDM8

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PRDM8 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PRDM8 antibody can be used for detection of PRDM8 by ELISA at 1:312500. PRDM8 antibody can be used for detection of PRDM8 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

72 kDa

Uniprot ID

[Q9NQV8](#)

NCBI

[NP_064611](#)

Dilution Range


PRDM8 antibody can be used for detection of PRDM8 by ELISA at 1:312500. PRDM8 antibody can be used for detection of PRDM8 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

Antibody used in WB on Human MCF-7 at 0....

PRDM8



Antibody used in WB on Human 293T at 1 u...