

Product Datasheet

XK Antibody (orb1246239)

Description	XK Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human XK.	Antibody used in WB on Human MCF-7 at 0....
Target	XK	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store XK antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	XK antibody can be used for detection of XK by ELISA at 1:12500. XK antibody can be used for detection of XK by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	51 kDa	
Uniprot ID	P51811	
NCBI	NP_066569	
Dilution Range	XK antibody can be used for detection of XK by ELISA at 1:12500. XK antibody can be used for detection of XK by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	