

Product Datasheet

TNKS Antibody (orb1246244)

Description

TNKS Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TNKS.

Target

TNKS

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TNKS antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TNKS antibody can be used for detection of TNKS by ELISA at 1:62500. TNKS antibody can be used for detection of TNKS by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

142 kDa

Uniprot ID

[O95271](#)

NCBI

[NP_003738](#)

Dilution Range

TNKS antibody can be used for detection of TNKS by ELISA at 1:62500. TNKS antibody can be used for detection of TNKS by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Small Inte...



Antibody used in WB on Human HepG2 at 1 ...



Antibody used in IHC on Human Small inte...