


Product Datasheet

KDM4C Antibody (orb1246296)

Description	KDM4C Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human JMJD2C.	Antibody used in WB on Human Jurkat 1 ug...
Target	KDM4C	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store JMJD2C antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	JMJD2C antibody can be used for detection of JMJD2C by ELISA at 1:62500. JMJD2C antibody can be used for detection of JMJD2C by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	120 kDa	
Uniprot ID	Q9H3R0	
NCBI	NP_055876	
Dilution Range	JMJD2C antibody can be used for detection of JMJD2C by ELISA at 1:62500. JMJD2C antibody can be used for detection of JMJD2C by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	