




Product Datasheet

MORF4L1 Antibody (orb1246306)

Description	MORF4L1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MORF4L1.	Antibody used in WB on Human Placenta at...
Target	MORF4L1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MORF4L1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Hum. Adult Placen...
Note	For research use only	
Application notes	MORF4L1 antibody can be used for detection of MORF4L1 by ELISA at 1:500. MORF4L1 antibody can be used for detection of MORF4L1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	37 kDa	
Uniprot ID	Q5NVP9	Antibody used in WB on Hum. Fetal Brain ...
NCBI	NP_006782	
Dilution Range	MORF4L1 antibody can be used for detection of MORF4L1 by ELISA at 1:500. MORF4L1 antibody can be used for detection of MORF4L1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	