

Product Datasheet

ATOH1 Antibody (orb1246341)


Description

ATOH1 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ATOH1.
Target	ATOH1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ATOH1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ATOH1 antibody can be used for detection of ATOH1 by ELISA at 1:62500. ATOH1 antibody can be used for detection of ATOH1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	Q92858
NCBI	NP_005163
Dilution Range	ATOH1 antibody can be used for detection of ATOH1 by ELISA at 1:62500. ATOH1 antibody can be used for detection of ATOH1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Placenta at...



Antibody used in IHC on mouse cochlea.



Antibody used in IHC on Human Spleen.