

Product Datasheet

ZC3H7B Antibody (orb1246372)

Description

ZC3H7B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
721_B at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZC3H7B.

Target

ZC3H7B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZC3H7B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZC3H7B antibody can be used for detection of ZC3H7B by ELISA at 1:312500. ZC3H7B antibody can be used for detection of ZC3H7B by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

110 kDa

Uniprot ID

[Q9UGR2](#)

NCBI

[NP_060060](#)

Dilution Range

ZC3H7B antibody can be used for detection of ZC3H7B by ELISA at 1:312500. ZC3H7B antibody can be used for detection of ZC3H7B by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.