





Product Datasheet

AGAP2 Antibody (orb1246404)

Description	AGAP2 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-QASLDSIREAVINSQ	orb1246404 (0.3 ug/ml) staining of Human...
Target	AGAP2	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1246404 (3.75 ug/ml) staining of para...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 100kDa band observed in Human Brain (Amygdala and Cerebral Cortex) lysates (calculated MW of 90.5kDa according to NP_055585.1). Recommended concentration: 0.3-1ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Thyroid and Kidney. Recommended concentration: 3.75ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q99490	orb1246404 (3.75 ug/ml) staining of para...
NCBI	NP_055585.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 100kDa band observed in Human Brain (Amygdala and Cerebral Cortex) lysates (calculated MW of 90.5kDa according to NP_055585.1). Recommended concentration: 0.3-1ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Thyroid and Kidney. Recommended concentration: 3.75ug/ml.	
Expiration Date	12 months from date of receipt.	