

Product Datasheet

EXOSC9 Antibody (orb1246431)

Description

EXOSC9 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

The immunogen for this antibody is: C-RTQTTSKQEKAP

Target

EXOSC9

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx. 75kDa band observed in lysates of cell line Jurkat (calculated MW of 50.8kDa according to NP_001029366.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Alderuccio, J Exp Med. 1991 Apr 1;173(4):941-52; PMID: 2007859). Recommended concentration: 0.3-1ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q06265](#)

NCBI


[NP_005024.2](#), [NP_001029366.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx. 75kDa band observed in lysates of cell line Jurkat (calculated MW of 50.8kDa according to NP_001029366.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Alderuccio, J Exp Med. 1991 Apr 1;173(4):941-52; PMID: 2007859). Recommended concentration: 0.3-1ug/ml.

Expiration Date

12 months from date of receipt.



orb1246431
(0.3 ug/ml)
staining of
Jurka...