

Product Datasheet

RASSF3 Antibody (orb1246502)

Description	RASSF3 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1246502 (0.5 ug/ml) staining of Heart...
Immunogen	The immunogen for this antibody is: CKPPQTSPNSGK	
Target	RASSF3	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:2000. Western Blot: Approx 26kDa band observed in Human Heart lysates (calculated MW of 27.6kDa according to NP_835463.1). Recommended concentration: 0.5-1.5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q86WH2	
NCBI	NP_835463.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:2000. Western Blot: Approx 26kDa band observed in Human Heart lysates (calculated MW of 27.6kDa according to NP_835463.1). Recommended concentration: 0.5-1.5ug/ml.	
Expiration Date	12 months from date of receipt.	