

## Product Datasheet

### NGFRAP1 Antibody (orb1246517)

## Description

NGFRAP1 Antibody

|                            |                                                                                                                                                                                                                                                                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Goat                                                                                                                                                                                                                                                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                                                                                                     |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen</b>           | The immunogen for this antibody is:<br>CLRILMGELSNHHDH                                                                                                                                                                                                                                                                                                           |
| <b>Target</b>              | NGFRAP1                                                                                                                                                                                                                                                                                                                                                          |
| <b>Preservatives</b>       | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                                      |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                                                                                                                                                                           |
| <b>Concentration</b>       | 500 ug/mL                                                                                                                                                                                                                                                                                                                                                        |
| <b>Storage</b>             | Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                                                                                                                                       |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                            |
| <b>Application notes</b>   | Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Preliminary experiments in Human Liver, Ovary, Prostate and Testis lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q00994</a>                                                                                                                                                                                                                                                                                                                                           |
| <b>NCBI</b>                | <a href="#">NP_996798.1</a> , <a href="#">NP_996800.1</a> , <a href="#">NP_055195.1</a>                                                                                                                                                                                                                                                                          |
| <b>Dilution Range</b>      | Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Preliminary experiments in Human Liver, Ovary, Prostate and Testis lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                                  |