


Product Datasheet

MTHFS Antibody (orb1246522)

Description		MTHFS Antibody		
Species/Host	Goat			
Reactivity	Human			
Conjugation	Unconjugated			
Tested Applications	ELISA, WB			
Immunogen	The immunogen for this antibody is: KRCLQHGEVKPYT		orb1246522 (0.5 ug/ml) staining of Human...	
Target	MTHFS			
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.			
Form/Appearance	Liquid			
Concentration	500 ug/mL			
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.			
Note	For research use only			
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 25kDa band observed in Human Liver lysates (calculated MW of 23.3kDa according to NP_006432.1). Recommended concentration: 0.5-1.5ug/ml.			
Clonality	Polyclonal			
Uniprot ID	P49914			
NCBI	NP_006432.1			
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 25kDa band observed in Human Liver lysates (calculated MW of 23.3kDa according to NP_006432.1). Recommended concentration: 0.5-1.5ug/ml.			
Expiration Date	12 months from date of receipt.			