


Product Datasheet

ABCC11 Antibody (orb1246525)

Description	ABCC11 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-KVVEFDRPEVLRK	
Target	ABCC11	orb1246525 (0.5 ug/ml) staining of Human...
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	orb1246525 (3.75 ug/ml) staining of para...
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx. 125kDa band observed in Human Breast and Human Testis lysates (calculated MW of 150kDa according to NP_660187.1). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.5-1.5ug/ml. Primary incubation 1 hour at room temperature. Immunohistochemistry: Paraffin embedded Human Breast. Recommended concentration: 3.75ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q96J66	
NCBI	NP_660187.1 , NP_115972.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx. 125kDa band observed in Human Breast and Human Testis lysates (calculated MW of 150kDa according to NP_660187.1). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.5-1.5ug/ml. Primary incubation 1 hour at room temperature. Immunohistochemistry: Paraffin embedded Human Breast. Recommended concentration: 3.75ug/ml.	
Expiration Date	12 months from date of receipt.	