

Product Datasheet

ETV1 Antibody (orb1246552)

Description

ETV1 Antibody

Species/Host

Goat

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

The immunogen for this antibody is:
PLKIKKEPHSPCSE

Target

ETV1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Preliminary experiments in lysates in cell line A431 gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Clonality

Polyclonal

Uniprot ID
P50549
NCBI

**NP_001156622.1, NP_001156623.1,
NP_001156620.1, NP_001156619.1,
NP_001156624.1, NP_004947.2,
NP_001156621.1**

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Preliminary experiments in lysates in cell line A431 gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0) 1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558