

Product Datasheet

ERAP2 Antibody [Out of stock] (orb1246554)

Description

ERAP2 Antibody [Out of stock]

250kDa
150kDa
100kDa

Species/Host

Goat

Conjugation

Unconjugated

Tested

Applications

ELISA, WB

50kDa
37kDa
25kDa
20kDa
15kDa

Immunogen

The immunogen for this antibody is: C-QYLKKFSYRNAKND

orb1246554
(0.5 ug/ml)
staining of
Splee...

Target

ERAP2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx. 110kDa band observed in Human Spleen lysates (calculated MW of 110kDa according to NP_071745.1). Recommended concentration: 0.5-1.5ug/ml. An additional band of unknown identity was also consistently observed at 50kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any such results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Clonality

Polyclonal

Uniprot ID

Q6P179

NCBI

NP_071745.1

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx. 110kDa band observed in Human Spleen lysates (calculated MW of 110kDa according to NP_071745.1). Recommended concentration: 0.5-1.5ug/ml. An additional band of unknown identity was also consistently observed at 50kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any such results been reported with other