

## Product Datasheet

### Sh2d1b/Sh2d1b2 (mouse) Antibody (orb1246555)

## Description

Sh2d1b/Sh2d1b2 (mouse) Antibody

## Species/Host

Goat

## Conjugation

Unconjugated

## Tested

ELISA

## Applications

## Immunogen

The immunogen for this antibody is: C-ELNVYENTDEEYVD

## Target

Sh2d1b / Sh2d1b2 (mouse)

## Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

## Form/Appearance

Liquid

## Concentration

500 ug/mL

## Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

## Note

For research use only

## Application notes

Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

## Clonality

Polyclonal

## Uniprot ID

[O35324](#)

## NCBI

[NP\\_036139.2](#)

## Dilution Range

Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been