

Product Datasheet

DMP1 Antibody (orb1246559)

Description		
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-ENSEQANTQDSG	orb1246559 (1 ug/ml) staining of Human K...
Target	DMP1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:4000.Western Blot:Approx. 55kDa band observed in Human Ileum and Kidney lysates (calculated MW of 55.8kDa according to NP_004398.1 and of 54.1 according to NP_001073380.1). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q13316	
NCBI	NP_001073380.1, NP_004398.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000.Western Blot:Approx. 55kDa band observed in Human Ileum and Kidney lysates (calculated MW of 55.8kDa according to NP_004398.1 and of 54.1 according to NP_001073380.1). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	