


Product Datasheet

DLG1 Antibody (orb1246560)

Description	DLG1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-NSKTRDKGEIPDD	orb1246560 (0.3 ug/ml) staining of KELLY...
Target	DLG1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 150kDa band observed in lysates of cell line KELLY (calculated MW of 103kDa according to NP_004078.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q12959	
NCBI	NP_001091894.1 , NP_004078.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 150kDa band observed in lysates of cell line KELLY (calculated MW of 103kDa according to NP_004078.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	