



Product Datasheet

ABCD3 Antibody (orb1246591)

Description	ABCD3 Antibody	
Species/Host	Goat	
Reactivity	Human	orb1246591 (0.2 ug/ml) staining of Human...
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-PDGREDQKRKGISD	
Target	ABCD3	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1246591 (3.8 ug/ml) staining of paraf...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx. 70kDa band observed in Human Kidney lysates (calculated MW of 75.5kDa according to NP_002849.1). Recommended concentration: 0.2-0.6 ug/ml. Primary incubation was 1 hour. Immunohistochemistry: In paraffin embedded Human Testis shows vesicular staining in Sertoli cells. Recommended concentration: 2.5-3.8ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P28288	
NCBI	NP_002849.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx. 70kDa band observed in Human Kidney lysates (calculated MW of 75.5kDa according to NP_002849.1). Recommended concentration: 0.2-0.6 ug/ml. Primary incubation was 1 hour. Immunohistochemistry: In paraffin embedded Human Testis shows vesicular staining in Sertoli cells. Recommended concentration: 2.5-3.8ug/ml.	
Expiration Date	12 months from date of receipt.	