




Product Datasheet

CSRP3 Antibody (orb1246610)

Description	CSRP3 Antibody	
Species/Host	Goat	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1246610 (0.01 ug/ml) staining of Huma...
Immunogen	The immunogen for this antibody is: CGKSLESTNVTDKD	
Target	CSRP3	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1246610 (2.5 ug/ml) staining of paraf...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx. 20kDa band observed in Human, Mouse and Rat Heart lysates (calculated MW of 21.0kDa according to Human NP_003467.1 and 20.8kDa according to Rat NP_476485.1). Recommended concentration: 0.01-0.03ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Skeletal Muscle and Heart. Recommended concentration: 2.5ug/ml.	
Clonality	Polyclonal	orb1246610 (2.5 ug/ml) staining of paraf...
Uniprot ID	P50461	
NCBI	NP_003467.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx. 20kDa band observed in Human, Mouse and Rat Heart lysates (calculated MW of 21.0kDa according to Human NP_003467.1 and 20.8kDa according to Rat NP_476485.1). Recommended concentration: 0.01-0.03ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Skeletal Muscle and Heart. Recommended concentration: 2.5ug/ml.	
Expiration Date	12 months from date of receipt.	