

Product Datasheet

UCN3 Antibody (orb1246662)

Description

UCN3 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

The immunogen for this antibody is: C-KPRQDTAKSPHRTK

Target

UCN3

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: In transfected HEK293 transiently expressing full-length Human UCN3 (myc and DYKDDDDK tagged), a band of approx. 19kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-DYKDDDDK tag antibody. Recommended concentration, 1- 2ug/ml.

Clonality


Polyclonal

Uniprot ID
[Q969E3](#)
NCBI
[NP_444277.2](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: In transfected HEK293 transiently expressing full-length Human UCN3 (myc and DYKDDDDK tagged), a band of approx. 19kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-DYKDDDDK tag antibody. Recommended concentration, 1- 2ug/ml.

Expiration Date

12 months from date of receipt.



HEK293
lysate (10 ug
protein in
RIPA buf...