

Product Datasheet

RPGRIP1L Antibody (orb1246674)

Description

RPGRIP1L Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

The immunogen for this antibody is: C-DKENNKAKRDILK

Target

RPGRIP1L

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 150kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 151kDa according to NP_056087.2). Recommended concentration: 0.3-1ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q68CZ1](#)

NCBI


[NP_056087.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 150kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 151kDa according to NP_056087.2). Recommended concentration: 0.3-1ug/ml.

Expiration Date

12 months from date of receipt.



orb1246674
(0.3 ug/ml)
staining of
Human...