

Product Datasheet

ZCCHC11 Antibody (orb1246769)

Description

ZCCHC11 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications

Immunogen

The immunogen for this antibody is: RERCPHPPRGNVSE

Target

ZCCHC11

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:16000. Immunohistochemistry: In paraffin embedded Human Lymph Node shows nuclear staining of select lymphocytes and staining of macrophages. Recommended concentration: 2-6ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q5TAX3](#)

NCBI


[NP_001009881.1](#), [NP_056084.1](#), [NP_001009882.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:16000. Immunohistochemistry: In paraffin embedded Human Lymph Node shows nuclear staining of select lymphocytes and staining of macrophages. Recommended concentration: 2-6ug/ml.

Expiration Date

12 months from date of receipt.



orb1246769
(4 ug/ml)
staining of
paraffi...