

Product Datasheet

LDB3 Antibody (orb1246770)

Description

LDB3 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1246770
(0.1 ug/ml)
staining of
Human...

Immunogen

The immunogen for this antibody is: C-KSQNKPEDEADE

Target

LDB3

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Western Blot experiments gave bands at approx 90kDa and 35kDa in Human Skeletal Muscle lysate after 1ug/ml antibody staining. These bands correspond to earlier findings in literature with different antibodies (Huang et al J Biol Chem. 2003 Feb 28;278(9):7360-5; PMID: 12499364). This protein has a calculated MW of 77.1kDa according to NP_009009.1. Recommended concentration: 0.1-1ug/ml.



orb1246770
(3.8 ug/ml)
staining of
paraf...

Clonality

Polyclonal

Uniprot ID

[O75112](#)

NCBI

[NP_001073584.1](#), [NP_001073585.1](#), [NP_009009.1](#),
[NP_001073583.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Western Blot experiments gave bands at approx 90kDa and 35kDa in Human Skeletal Muscle lysate after 1ug/ml antibody staining. These bands correspond to earlier findings in literature with different antibodies (Huang et al J Biol Chem. 2003 Feb 28;278(9):7360-5; PMID: 12499364). This protein has a calculated MW of 77.1kDa according to NP_009009.1. Recommended concentration: 0.1-1ug/ml.

Expiration Date

12 months from date of receipt.