




Product Datasheet

VPS28 Antibody (orb1246793)

Description	VPS28 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1246793 staining (1 ug/ml) of Human T...
Immunogen	The immunogen for this antibody is: C-ESAYNAFNRLHA	
Target	VPS28	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1246793 (10 ug/ml) staining of paraff...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Approx 24kDa band observed in H460, Human Testis and Brain lysates (predicted MW of 25kDa according to NP_057292.1). Recommended for use at 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Testis, stains acrosomal head caps in round spermatids. Recommended concentration: 10ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9UK41	
NCBI	NP_057292.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Approx 24kDa band observed in H460, Human Testis and Brain lysates (predicted MW of 25kDa according to NP_057292.1). Recommended for use at 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Testis, stains acrosomal head caps in round spermatids. Recommended concentration: 10ug/ml.	
Expiration Date	12 months from date of receipt.	